
DE - Produktdetails

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EN - Product Details

- With its protected shape, the cushion especially scores as a dynamic bolster and back support cushion. The flowing transition between the cushion and the seat surface prevents press marks at the edge.
- In an erect position, your vertebral disks move gently. Minimal movements stimulate the deep muscles. Your back is supported, relieved and massaged.
- Strong abdominal muscles stabilise the spine and are especially important when it comes to preventing back complaints. During an abdominal workout, the flat and unstable cushion can be unobtrusively positioned under your buttocks. The deep muscles are additionally activated and your abdominal training is more intensive.
- Specially formed knobs permit an activation of your foot reflex zones.

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EN - Sitting Recommendations

- Modern sitting research recommends movement and variety when sitting. Introduce variety to your daily sitting marathon and occasionally use the cushion as a dynamic sitting support or as a back support cushion.
- When doing so, you should adapt the sitting duration or sitting frequency on the cushion to your fitness level.
- The closer you sit to the front of the cushion, the lower the sideways movement of the cushion. Such a sitting position is recommended if you want to prolong the sitting phases on the cushion without placing excessive burden on the back muscles.

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EN - Air Regulation

- Remove the needle from the pump handle on the side!
- Insert the needle into the valve as far as possible. Important: take hold of the valve with your fingers and twist it outwards. You can now insert the needle into the valve diagonally in the direction of the centre of the cushion, without damaging the other side of the cushion.
- To reduce the air quantity, press down on the cushion so that air can flow out via the needle (without the pump). Air is added using the pump.
- The valve should be positioned in the rear part of the cushion so that you do not notice it when there is only a small amount of air in the cushion. Turn the cushion so that it is correctly positioned before you sit on it!

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More Exercises

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Foot Reflexology